



# Specification Information

**Product Code:** 878-500

**Product Name:** Internal Disposable Containers for Code 878

**Individual Carton UPC:** 624622878054

**Individual Carton Shipping Dimensions:** 50.8cm L x 17.8cm W x 12.7cm H (20" L x 7" W x 5" H)

**Individual Carton Shipping Weight:** 1.8kg (4 lbs)

**Site Location:** Frost's 878-500 internal disposable containers are designed for use in code 878. These units are perfect for health care and commercial environments. Typically, this product is used in hospitals, long term care facilities, doctor/dentist offices, schools/universities, office buildings, and retail locations. The high quality and rigorous testing of this product make it a popular choice.



## Features:

- Case of four internal plastic disposals per box
- Compatible with Frost code 878 providing optimal convenience
- Suitable for discarding of: needles, syringes, blades, pins, gauze, and cotton balls
- Internal containers are designed and certified to exceed penetration requirements for all applicable international standards (FDA #9613648)

## Capacity:

- Each plastic disposable container has a capacity of 4 litres (4.23qt)

## Materials:

- High quality plastic

## Product Dimensions:

- 14.0cm W x 30.2cm H x 14.0cm D (5.5" W x 11.9" H x 5.5" D)

## Product Weight:

- 0.13 kg. (0.3 lbs) each, case weighs 0.54 kg. (1.2 lbs)

**Installation Guidelines:** For use with Frost code 878, place internal plastic container's mouth into the bottom of the inlet tube and slide the base of the container onto the bottom of the holder. Please ensure the lid is to either side of the plastic disposable container to ensure the door closes properly. Please note that it is recommended that all installations

Frost Products Ltd

5276 John Lucas Drive  
Burlington, ON, Canada, L7L5Z9

(905) 332-5700  
1-800 667-5094

info@frostproductsltd.com  
www.frostproductsltd.com



## *Specification Information*

be conducted by a qualified tradesperson. Please read and follow any additional installation guidelines provided with the products.

**Maintenance:** Maintenance staff should periodically inspect internal container's fullness based on site usage. If it has reached capacity, then remove and dispose of according to local guidelines. Replace with unused internal container (code 878-500).

**Warranty:** If any product or part proves to be defective within 12 months from date of shipment from the factory it will be repaired or replaced without any charge including freight costs when authorized. Costs resulting from the installation or removal of defective products are not covered under this warranty. Discolouration to stainless steel parts resulting from exposure to harsh environmental and / or chemicals are not considered to be defects in workmanship or material and are not covered under warranty.